
**Agilent 3200SA
Stirrer
搅拌器**

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Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

5973-1774

Edition

First Edition, March 2012

Printed in China

Agilent Technologies, Inc.
412 Ying Lun Road
Shanghai 200131 P.R. China

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Safety Notices

CAUTION


A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1	Installation	5
	Tools and Accessories	6
	Installing the Stirrer	7
	Installing the electrode pole	7
	Installing the electrode fixture	8
	Installing the power adapter	9
2	Operation	11
	Introduction	12
	Where to find information	12
	Agilent customer portal	12
	Features	13
	Specifications	13
	Uses for the 3200SA Stirrer	14
	Using the 3200SA Stirrer	15
	Prepare a container	15
	Install electrodes	16
	Start the 3200SA Stirrer	16
	Turn off the 3200SA Stirrer	16
	Operation panel	17
	Display	17
	Keypad	17
3	Troubleshooting and Maintenance	19
	Troubleshooting	20
	Maintenance	20
	Consumables	21

1 Installation

Tools and Accessories	6
Installing the Stirrer	7

Tools and Accessories

Agilent provides the tools needed for installation. The following parts are included in the shipping case.

- Electrode pole
- B620 Stirrer bar (G4388-27000)
- Electrode fixture
- Power adapter (5185-8389)

Installing the Stirrer

Open the 3200SA Stirrer shipping case. Remove the stirrer and other accessories.

Installing the electrode pole

Thread the electrode pole into the top of the stirrer as shown in [Figure 1](#).

Figure 1 Installing the electrode pole

Installing the electrode fixture

- 1 Hold the button on the side of the electrode fixture and slide the fixture onto the electrode pole.
- 2 Release the button when the fixture is in place. See [Figure 2](#).


Figure 2 Installing the electrode fixture

Installing the power adapter

A universal power adapter is included with the 3200SA. The power adapter operates at 100–240 VAC, 1 A, 50/60 Hz. The output from the power adapter is 9 VDC 1 A.

CAUTION

Do not use any other power adapter with the 3200SA.

- 1 The power adapter includes several outlet adapters. Choose the appropriate outlet adapter for the power outlets in your country.
- 2 Slide the outlet adapter on to the two metal prongs on the power adapter. A click will be heard when the plug is properly engaged. See [Figure 3](#).


Figure 3 Assembling the power adapter

1 Installation

- 3 Connect the power adapter cord to the power socket on the back of the stirrer. See [Figure 4](#).

Figure 4 Back view of the stirrer

2 Operation

Introduction	12
Using the 3200SA Stirrer	15

Introduction

Where to find information

This document describes how to install, operate, maintain, and troubleshoot the 3200SA Stirrer.

Read this document before operating the stirrer. To avoid common safety issues:

- Set up the stirrer in a spacious area on an even, stable, clean, non-slip, dry, and fireproof surface.
- Check the instrument before use, do not use a damaged instrument.
- Keep the base plate clean.

CAUTION

The stirrer may heat up when in use. Do not cover the stirrer, even partially, to avoid overheating.

CAUTION

Only use the power adapter included with the stirrer.

WARNING

Wear personal protective equipment. Stirring may splash or vaporize liquids, propel an object out of the liquid container, or release toxic or combustible gases. Hazards may occur from flammable materials, combustible liquid with a low boiling temperature, incorrect liquid container size, or general unsafe container conditions.

Agilent customer portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at <http://www.agilent.com/chem>.

Features

The 3200SA Stirrer consists of a base containing the controls and magnetic coil, a fixture to hold electrodes in place, and a stirring bar. It can be used with various measurement electrodes, and it can be used to stir any liquids.

The 3200SA Stirrer operates with an advanced magnetic coil. There are no moving parts inside the stirrer. The rotating magnetic stirring bar is driven by alternately changing the polarity of the internal magnetic coil. This works because similar magnetic poles repel each other.

The stirring bar can rotate at a wide range of speeds, and once the speed is set it rotates at a steady speed.

Specifications

Speed	0–1500 rpm
Stirred volume (max)	1000 mL (at 1500 rpm with 6 mm diameter, 20 mm long stirring bar)
Stirring bar	30 × 8 mm (max)
Environmental temperature	0 to 40 °C
Environmental relative humidity	≤ 85%
Power supply	Power adapter (5185-8389)
Power input	100–240 VAC, 1 A
Power output	9 VDC, 1 A
Dimensions (mm)	200 length × 120 width × 50 height
Weight	approximately 0.5 kg

CAUTION

Do not use where nearby vibrations will affect the performance of the stirrer. Do not use if corrosive gas is in the air.

Uses for the 3200SA Stirrer

The stirrer can be used in university laboratories and industrial enterprises to stir sample solutions. It can also be used with electrodes to make measurements.

Using the 3200SA Stirrer

CAUTION

The stirrer must be kept clean and dry. Avoid contact of the socket with acidic, alkaline, and salt solutions.

Prepare a container

Prepare a liquid container before using the stirrer. To prepare a container:

- 1 Place the stirring bar in a container. Magnetic fields must be able to penetrate the container.
- 2 Fill the container with the sample solution to be stirred.
- 3 Place the container in the center of the stirring base so the stirring bar will rotate smoothly. Use the circle engraved on the stirring base as reference. See [Figure 5](#).


Figure 5 Front view of the stirrer

Install electrodes

Install measurement electrodes if needed.

- 1 Attach the electrode to the electrode fixture. See [Figure 5 on page 15](#).
- 2 Hold the button on the side of the electrode fixture to adjust the height. Keep enough distance between the tip of the electrode and the stirring bar.

Start the 3200SA Stirrer

The 3200SA Stirrer is controlled by the keypad shown in [Figure 6 on page 17](#).

- 1 Press **[ON/OFF]** to turn on the stirrer and start stirring. The screen displays the stirring speed.
- 2 Press **▲** and **▼** keys to adjust the speed of the stirrer bar. The speed setting remains unchanged when the stirrer is unplugged.

NOTE

At certain speeds, irregularities in the shape of the stirring bar and the bottom of the liquid container can cause the bar to move irregularly during stirring. Increase or decrease the speed until the stirring bar rotates evenly.

If the stirring bar does not rotate, check whether the container bottom is flat or if it obstructs the stirring bar. The container may also be severely off the center of the base plate. Adjust the location of the container.

Turn off the 3200SA Stirrer

After using the stirrer, press **[On/Off]** to turn off the stirrer. If the stirrer will be out of use for a long time, disconnect the power adapter from the back of the stirrer.

Operation panel

The operation panel consists of the display and the keypad.
See [Figure 5 on page 15](#).

Display

The display shows the current speed of the 3200SA Stirrer.

Keypad

The meter has 3 operating keys.



Figure 6 The keypad

3 Troubleshooting and Maintenance

Troubleshooting	20
Maintenance	20

Troubleshooting

Display does not light up	Check that the power adapter is connected. Use a volt meter to check the power output of the power adapter.
Irregular bar movement	Increase or decrease the speed until the stirring bar rotates evenly.
Bar does not rotate	Check the bottom of the container. A raised bottom may obstruct the movement of the bar. Check that the container is centered in the circle on the stirring base and adjust as needed.

Maintenance

CAUTION

Disconnect the power supply before cleaning. Wipe the stirrer with a soft cloth during cleaning. Do not allow moisture to get into the stirrer when cleaning.

WARNING

Wear protective gloves during cleaning.

Polluting substance	Remove using
Dyes	Isopropyl alcohol
Construction materials	Water with detergent Isopropyl alcohol
Cosmetics	Water with detergent Isopropyl alcohol
Foodstuffs	Water with detergent
Fuels	Water with detergent


Consumables

Part no.	Model and name	Description
5185-8389	Power adapter	100–240 VAC, 1 A, 50/60 Hz
G4388-27000	B620 Stirring bar	For use with the stirrer

Agilent provides customized information for the products you own through a customer portal. Log onto the portal at <http://www.agilent.com/chem>.

For new purchases, contact distributors of Agilent Technologies or log onto the Agilent official website.

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